



## Product Specification Summary

### Gene Detection Systems

Catalog Number	40-2004-40
Product Name	Fragile X Syndrome GeneProber unlabeled Probe; 500 ng
Size	500 ng
Description	Fragile X unlabeled probe
Component/Note	CGG triple repeat spanning
Component/Note	Fragile X CGG
Storage	Store at -20°C

Fragile X GeneProber GLFX1 unlabeled Probe [Old Catalog Number 40-3201-01]. Fragile X CGG triple repeat spanning region unlabeled probe for radioactive labeling and Southern blot detection. Suitable for random primer labeling. 500ng lyophilized product For research use only. Not for use in diagnostic procedures for clinical purposes

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=74>



Product Manual

<http://www.genelink.com/Literature/ps/M40-2004-40Ver5.3.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Factor V Leiden Genotyping Kit; 1 Kit 100 Reactions	40-2035-10K	1 Kit
Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit	40-2036-11	1 Kit
Custom Genotyping	40-2000-01	1 Ea
ALS C9orf72 hexanucleotide expansion Genotyping GeneProber ALS-GL577	40-2057-41	110 uL
Factor VIII Genemer Pack (Hemophilia A); 10 nmols	40-2036-10	10 nmols
Factor V Leiden Genemer; 2 nmols	40-2035-10	2 nmols
X alphoid repeat Genemer; 10 nmoles	40-2021-10	10 nmols
Y alphoid repeat Genemer; 10 nmols	40-2022-10	10 nmols
Sickle Cell Genemer HbA control DNA; 500 ng	40-2001-01	500 ng
Sickle Cell Genemer HbS control DNA; 500 ng	40-2001-02	500 ng
Sickle Cell Genemer HbC control DNA; 500 ng	40-2001-03	500 ng
Fragile X CGG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL	40-2004-41	110 uL
Fragile X Syndrome GeneProber unlabeled Probe; 500 ng	40-2004-40	500 ng
Fragile X Syndrome PCRProber Kit; 5 blots	40-2004-32	5 blots
GLFX Genemer Kit for Radioactive Detection; 1 run	40-2004-20	1 run