



Product Specification Summary

GeneProber™ Digoxigenin Labeled Probes

| | |
|----------------|--|
| Catalog Number | 40-2025-41 |
| Product Name | Huntington CAG Repeat Genotyping GeneProber™ Digoxigenin labeled; 110 uL |
| Size | 110 uL |
| Description | Huntington CAG Repeat Genotyping GeneProber™ Digoxigenin labeled; 110 uL |
| Component/Note | CAG Triple repeat region |
| Component/Note | Huntington CAG Repeat Genotyping GeneProber™ Digoxigenin labeled; 110 uL |

Gene Link is the leader in triple repeat disorder genotyping using non-radioactive based methods. We have more than a decade of expertise and have developed facile non-radioactive detection methods for safe, sensitive and reliable genotyping of human genetic disorders specifically fragile X, Huntington Disease, Myotonic Dystrophy etc. Take a look at the simple agarose and polyacrylamide gel based systems, the chemiluminescent Southern blot detection methods and the fluorescent systems for genotyping of triple repeat disorders.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=864>



Product Manual

http://www.genelink.com/Literature/ps/M40-2025-41_Ver6.1.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|--------|
| Kennedy Disease (SBMA) GeneProber Digoxigenin labeled Probe; 100 uL | 40-2032-41 | 110 uL |
| SCA7 GeneProber Digoxigenin labeled Probe; 110 uL | 40-2041-41 | 110 uL |
| DRPLA GeneProber Digoxigenin labeled; 110 uL | 40-2042-41 | 110 uL |
| Human mtDNA Deletion Genotyping GeneProber GL557 Digoxigenin labeled ; 110 uL | 40-2055-41 | 110 uL |
| ALS C9orf72 hexanucleotide expansion Genotyping GeneProber ALS-GL577 | 40-2057-41 | 110 uL |
| Mouse 18S rDNA genomic GeneProber Digoxigenin Labeled | 40-2058-41 | 110 uL |
| SCA1 GeneProber Digoxigenin labeled Probe; 110 uL | 40-2037-41 | 110 uL |
| SCA2 GeneProber Digoxigenin labeled Probe; 110 uL | 40-2038-41 | 110 uL |
| SCA3 GeneProber Digoxigenin labeled Probe; 110 uL | 40-2039-41 | 110 uL |
| SCA6 GeneProber Digoxigenin labeled Probe; 110 uL | 40-2040-41 | 110 uL |
| Mouse mtDNA Genotyping GeneProber GL576 Digoxigenin Labeled | 40-2056-41 | 110 uL |
| Huntington CAG Repeat Genotyping GeneProber™ Digoxigenin labeled; 110 uL | 40-2025-41 | 110 uL |
| FRAXE/FMR2/AFF2 CCG Repeat Genotyping GeneProber™ AFF2-AJ31Dig1 | 40-2054-41 | 110 uL |
| Fragile X CGG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL | 40-2004-41 | 110 uL |
| Myotonic Dystrophy CTG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL | 40-2026-41 | 110 uL |