



Product Specification Summary

FMR2/AFF2/FRAXE

Catalog Number	40-2054-41
Product Name	FRAXE/FMR2/AFF2 CCG Repeat Genotyping GeneProber™ AFF2-AJ31Dig1
Size	110 uL
Description	FRAXE/FMR2/AFF2 CCG Repeat Genotyping GeneProber™ AFF2-AJ31Dig1
Component/Note	CCG Triple repeat region
Component/Note	FRAXE/FMR2/AFF2 CCG Repeat Genotyping GeneProber™ AFF2-AJ31Dig1
Storage	Store at -20°C

We have developed a probe to be used for Southern blot analysis that reliably detects the AFF2 CCG triple repeat amplification. The presence of pre-mutated and fully expanded alleles in either gender, were confirmed by Southern blot analysis, which also enabled exclusion of methylation or repeat number mosaics as well as PCR failure. The use of this probe has been recommended as suitable for genotyping of pre-mutations, full mutations, and mosaics specifically for individuals presented for FRAXA screening with negative results to determine FRAXE status.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17177>



Product Manual

http://www.genelink.com/Literature/ps/M40-2054-41_V2.1.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
FRAXE/FMR2/AFF2 GGscan Kit CCG Repeat Genotyping; 100 rxns.	40-2054-15FM	1 Kit
FRAXE/FMR2/AFF2 GGscan Kit CCG Repeat Genotyping. 20 rxns	40-2054-15FMS	1 kit
FRAXE/FMR2/AFF2 CCG Repeat Genotyping GeneProber™ AFF2-AJ31Dig1	40-2054-41	110 uL