



Product Specification Summary

mtDNA Deletion

Catalog Number 40-2056-41
Product Name Mouse mtDNA Genotyping GeneProber GL576 Digoxigenin Labeled
Size 110 uL
Description Mouse Mitochondrial DNA ~3 kb span GeneProber GL576 Digoxigenin Labeled Probe
Component/Note **mtDNA Deletion GeneProber**
The probe spans ~ 3kb of mouse mt-DNA.
Component/Note One tube containing 110 µL of mouse mitochondrial DNA GeneProber GL576 digoxigenin labeled probe at a concentration of ~40ng/uL. The quantity supplied is sufficient for at least 5 20x20 cm blots using 20µL for each blot a
Scan the QR Code or visit the following links
Product Infomation <http://www.genelink.com/geneprodsite/product.asp?p=17359>



Product Manual http://www.genelink.com/Literature/ps/M40-2055-41_Ver1.3.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Human mtDNA Deletion Genotyping GeneProber GL557 Digoxigenin labeled ; 110 uL	40-2055-41	110 uL
Mouse mtDNA Genotyping GeneProber GL576 Digoxigenin Labeled	40-2056-41	110 uL