



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-3000-05
Product Name	DEPC Treated Water; 50 mL
Size	50 mL
Description	DEPC Treated Water; 50 mL

Gene Link DEPC treated water and Nuclease Free water (DEPC free) has been treated and autoclaved water that are certified to be free of nucleases and nucleic acids. Nuclease Free water (DEPC free) has not been treated with DEPC and is certified for all molecular biology applications requiring nuclease and nucleic acid free water. Quality Control includes RT-PCR, standard PCR and RNA gels without any degradation certifies the product to free of non-specific endonuclease, exonuclease and RNase activity.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17182>



Product Manual [http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA\\_Soln\\_Ver2.1.pdf](http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA_Soln_Ver2.1.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_40-DNA&RNA\\_solutions.pdf](http://www.genelink.com/Literature/ps/MSDS_40-DNA&RNA_solutions.pdf)



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL