

# Product Specification Summary

## DNA & RNA Reconstitution Solutions

Catalog Number 40-3000-05  
Product Name DEPC Treated Water; 50 mL  
Size 50 mL  
Description DEPC Treated Water; 50 mL

Gene Link DEPC treated water and Nuclease Free water (DEPC free) has been treated and autoclaved water that are certified to be free of nucleases and nucleic acids. Nuclease Free water (DEPC free) has not been treated with DEPC and is certified for all molecular biology applications requiring nuclease and nucleic acid free water. Quality Control includes RT-PCR, standard PCR and RNA gels without any degradation certifies the product to free of non-specific endonuclease, exonuclease and RNase activity.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17182>



Product Manual [http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA\\_Soln\\_Ver2.1.pdf](http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA_Soln_Ver2.1.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_40-DNA&RNA\\_solutions.pdf](http://www.genelink.com/Literature/ps/MSDS_40-DNA&RNA_solutions.pdf)



### Related Products

Product	Catalog No	Size
Nuclease Free Water (DEPC Free) 10 X 1.6 mL; 16 mL	40-3001-16	16 mL
DEPC Treated Water; 10 X 1.6 mL	40-3000-16	1.6 mL
DEPC Treated Water; 50 mL	40-3000-05	50 mL