



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-3000-16
Product Name	DEPC Treated Water; 10 X 1.6 mL
Size	1.6 mL
Description	DEPC Treated Water; 10 X 1.6 mL

Gene Link DEPC treated water and Nuclease Free water (DEPC free) has been treated and autoclaved water that are certified to be free of nucleases and nucleic acids. Nuclease Free water (DEPC free) has not been treated with DEPC and is certified for all molecular biology applications requiring nuclease and nucleic acid free water. The smaller 1.6 mL product size is specifically recommended for RT-PCR (QPCR) applications. Quality Control includes RT-PCR, standard PCR and RNA gels without any degradation certifies the product to free of non-specific endonuclease, exonuclease and RNase activity.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17181>



Product Manual [http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA\\_Soln\\_Ver2.1.pdf](http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA_Soln_Ver2.1.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_40-DNA&RNA\\_solutions.pdf](http://www.genelink.com/Literature/ps/MSDS_40-DNA&RNA_solutions.pdf)



### Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit

Omni-Marker™ Universal unlabeled; 100 uL

40-3005-01

190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |

SS40-3000-16.pdf

Print Date Version : May 23, 2017

Page 1 of 1

