



Product Specification Summary

Genetic Tools & Reagents

Catalog Number	40-3004-10
Product Name	Loading Buffer 5X Orange G/XC non-denaturing; 1 mL
Size	1 mL
Description	Orange G and XC dyes
Component/Note	NON-DENATURING
Storage	Store at 15-25°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=921>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3004-XXver1.2.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL