



# Product Specification Summary

## Genetic Tools & Reagents

|                |  |
|----------------|--|
| Catalog Number | 40-3004-10   |
| Product Name   | Loading Buffer 5X Orange G/XC non-denaturing; 1 mL |
| Size           | 1 mL   |
| Description    | Orange G and XC dyes                               |
| Component/Note | NON-DENATURING                                     |
| Storage        | Store at 15-25°C                                   |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=921>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3004-XXver1.2.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product  | Catalog No | Size   |
|--|------------|--------|
| Omni-Array™ mRN Antisense strand amplification kit | 08-0021-02 | 1 Kit  |
| Omni-Array™ Sense Strand mRNA Amplification Kit    | 08-0011-02 | 1 Kit  |
| Omni-Marker™ Universal unlabeled; 100 uL           | 40-3005-01 | 100 uL |