



Product Specification Summary

Omni-Marker

Catalog Number	40-3005-01
Product Name	Omni-Marker™ Universal unlabeled; 100 uL
Size	100 uL
Description	50 bp to 10 kb Fragments
Component/Note	Load 5 to 10 ul
Storage	Store at -20°C

Omni-Marker™ Universal and Low are unlabeled DNA markers containing a blend of fragments ranging from 50 base pairs to 10 kb. The universal contains fragments of the following sizes; 10 kb, 8 kb, 6 kb, 4 kb, 3 kb, 2 kb, 1.55, 1.4 kb, 1 kb, 750 bp, 500 bp, 400 bp, all the bands and "low" version contains fragments from 50 bp to 2kb. The low Omni-Marker™ is ideal for routine PCR gels. These are ready to load containing BPB/XC dyes. A loading of 5 µl is sufficient per lane.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=900>



Product Manual

<http://www.genelink.com/Literature/ps/Omni-Marker ver2.3.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni- Marker™ Universal unlabeled; 500 uL	40-3005-05	500 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni- Marker™ Low unlabeled; 100 uL	40-3006-01	100 uL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL
Omni- Marker™ Low unlabeled; 1 mL	40-3006-10	1 mL
GScan Marker Tamra Labeled 50 bp - 600 bp; 100 uL	40-3061-01	100 uL
GScan Marker Tamra Labeled 50 bp - 600 bp; 500 uL	40-3061-05	500 uL
Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb	40-3071-01	100 ul
Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb	40-3071-05	500 ul
Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp	40-3072-01	100 ul
Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp	40-3072-05	500 ul
GScan Marker Hex Labeled 50 bp - 600 bp; 100 uL	40-3081-01	100 uL
GScan Marker Hex Labeled 50 bp - 600 bp; 500 uL	40-3081-05	500 uL
Omni-Marker™ Digoxigenin labeled 50 bp-1kb	40-3091-01	100 ul