



Product Specification Summary

Genetic Tools & Reagents

| | |
|----------------|---|
| Catalog Number | 40-3005-05 |
| Product Name | Omni- Marker™ Universal unlabeled; 500 uL |
| Size | 500 uL |
| Description | 50 bp to 10 kb Fragments |
| Component/Note | Load 5 to 10 ul |
| Storage | Store at -20°C |

Omni-Marker™ Universal and Low are unlabeled DNA markers containing a blend of fragments ranging from 50 base pairs to 10 kb. The universal contains fragments of the following sizes; 10 kb, 8 kb, 6 kb, 4 kb, 3 kb, 2 kb, 1.55, 1.4 kb, 1 kb, 750 bp, 500 bp, 400 bp, all the bands and "low" version contains fragments from 50 bp to 2kb. The low Omni-Marker™ is ideal for routine PCR gels. These are ready to load containing BPB/XC dyes. A loading of 5 µl is sufficient per lane.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=901>



Product Manual

<http://www.genelink.com/Literature/ps/Omni-Marker ver2.3.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|--------|
| Omni-Array™ mRN Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |
| Omni- Marker™ Universal unlabeled; 500 uL | 40-3005-05 | 500 uL |
| Omni-Marker™ Universal unlabeled; 1 mL | 40-3005-10 | 1 mL |
| Omni- Marker™ Low unlabeled; 100 uL | 40-3006-01 | 100 uL |
| Omni-Marker™ Low unlabeled; 500 uL | 40-3006-05 | 500 uL |
| Omni- Marker™ Low unlabeled; 1 mL | 40-3006-10 | 1 mL |
| GScan Marker Tamra Labeled 50 bp - 600 bp; 100 uL | 40-3061-01 | 100 uL |
| GScan Marker Tamra Labeled 50 bp - 600 bp; 500 uL | 40-3061-05 | 500 uL |
| Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb | 40-3071-01 | 100 ul |
| Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb | 40-3071-05 | 500 ul |
| Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp | 40-3072-01 | 100 ul |
| Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp | 40-3072-05 | 500 ul |
| GScan Marker Hex Labeled 50 bp - 600 bp; 100 uL | 40-3081-01 | 100 uL |