



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-3005-10
Product Name	Omni-Marker™ Universal unlabeled; 1 mL
Size	1 mL
Component/Note	Load 5 to 10 ul

Omni-Marker™ Universal and Low are unlabeled DNA markers containing a blend of fragments ranging from 50 base pairs to 10 kb. The universal contains fragments of the following sizes; 10 kb, 8 kb, 6 kb, 4 kb, 3 kb, 2 kb, 1.55, 1.4 kb, 1 kb, 750 bp, 500 bp, 400 bp, all the bands and "low" version contains fragments from 50 bp to 2kb. The low Omni-Marker™ is ideal for routine PCR gels. These are ready to load containing BPB/XC dyes. A loading of 5 µl is sufficient per lane.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=902>



Product Manual <http://www.genelink.com/Literature/ps/Omni-Marker ver2.3.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL