



Product Specification Summary

PCR Amplification & Analysis

| | |
|----------------|--------------------------------------|
| Catalog Number | 40-3007-01 |
| Product Name | 50 X TAE buffer; 100 mL |
| Size | 100 mL |
| Description | TAE Buffer; 50 X Concentrate; 100 ml |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1072>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3007-XX Ver 1.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|--------|
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |
| Omni-Marker™ Universal unlabeled; 1 mL | 40-3005-10 | 1 mL |
| Omni-Marker™ Low unlabeled; 500 uL | 40-3006-05 | 500 uL |