



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-3007-10
Product Name	TAE Buffer 50 X Concentrate; 1000 mL
Size	1000 mL
Description	TAE Buffer; 50 X Concentrate; 1000 ml

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1073>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3007-XX Ver 1.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL