



# Product Specification Summary

## PCR Amplification & Analysis

|                |                                       |
|----------------|---------------------------------------|
| Catalog Number | 40-3007-10                            |
| Product Name   | TAE Buffer 50 X Concentrate; 1000 mL  |
| Size           | 1000 mL                               |
| Description    | TAE Buffer; 50 X Concentrate; 1000 ml |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1073>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3007-XX Ver 1.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product                                  | Catalog No | Size   |
|--|------------|--------|
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |
| Omni-Marker™ Universal unlabeled; 1 mL   | 40-3005-10 | 1 mL   |
| Omni-Marker™ Low unlabeled; 500 uL       | 40-3006-05 | 500 uL |