



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-3010-01
Product Name	Agarose LE Molecular Biology Grade 1 Kg; 1000 g
Size	1000 g
Description	Genetic Tools & Reagents
Component/Note	Store in a tightly closed container
Component/Note	Certified Ultra Pure molecular biology grade DNase and Rnase-free
Storage	Store at Room Temperature

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=16961>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3010-XX Ver 1.2.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	