



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-3023-16
Product Name	MgCl <sub>2</sub> 50 mM; 1.6 mL
Size	1.6 mL
Description	Magnesium Chloride 50 mM
Component/Note	For PCR Optimization
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=11668>



Product Manual [http://www.genelink.com/Literature/ps/M40-3021-PCR\\_Additives\\_Ver5.1.pdf](http://www.genelink.com/Literature/ps/M40-3021-PCR_Additives_Ver5.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL