

Product Specification Summary

PCR Additives

Catalog Number	40-3023-16
Product Name	MgCl ₂ 50 mM; 1.6 mL
Size	1.6 mL
Description	Magnesium Chloride 50 mM
Component/Note	For PCR Optimization
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=11668>



Product Manual http://www.genelink.com/Literature/ps/M40-3021-PCR_Additives_Ver5.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
PCR Master Mix (2X) 100 reactions (2 x 1.3 ml); 100 rxns	40-5213-01	100 rxns
PCR Master Mix (2X) 200 reactions (4 x 1.3 ml); 200 rxns	40-5213-02	200 rxns
PCR Buffer Standard (10X); 1.6 mL	40-3060-16	1.6 mL