



Product Specification Summary

PCR Amplification & Analysis

Catalog Number	40-3059-10
Product Name	KCl 300 mM; 1 mL
Size	1 mL
Description	PCR Additive
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1068>



Product Manual http://www.genelink.com/Literature/ps/M40-3021-PCR_Additives_Ver5.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL