



Product Specification Summary

PCR Amplification & Analysis

Catalog Number	40-3061-16
Product Name	PCR Buffer Mg Free (10 X); 1.6 mL
Size	1.6 mL
Description	PCR Buffer without Mg ⁺⁺
Component/Note	of Mg ⁺⁺
Component/Note	Add desired concentration
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1060>



Product Manual http://www.genelink.com/Literature/ps/PS40-3060-XX_Ver3.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	