



Product Specification Summary

PCR Amplification & Analysis

| | |
|----------------|-------------------------------------|
| Catalog Number | 40-3061-16 |
| Product Name | PCR Buffer Mg Free (10 X); 1.6 mL |
| Size | 1.6 mL |
| Description | PCR Buffer without Mg ⁺⁺ |
| Component/Note | of Mg ⁺⁺ |
| Component/Note | Add desired concentration |
| Storage | Store at -20°C |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1060>



Product Manual http://www.genelink.com/Literature/ps/PS40-3060-XX_Ver3.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|--------|
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |
| Omni-Marker™ Universal unlabeled; 1 mL | 40-3005-10 | 1 mL |
| Omni-Marker™ Low unlabeled; 500 uL | 40-3006-05 | |