



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-3070-10
Product Name	Taq Polymerase Dilution Buffer; 1 mL
Size	1 mL
Description	Storage/ Dilution Buffer
Component/Note	TaqDNA Pol Storage Buffer
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1061>



Product Manual [http://www.genelink.com/Literature/ps/PS40-3060-XX\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/PS40-3060-XX_Ver3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL