



Product Specification Summary

Omni-Pure™ Genomic DNA Purification Systems

Catalog Number	40-4011-10
Product Name	GD1 Cell Suspension Solution 1000 ml
Size	1000 mL
Description	Omni-Pure DNA Purification System

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=11609>



Product Manual http://www.genelink.com/Literature/ps/GDComponents_40-4010-XX.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf



Related Products

Product	Catalog No	Size
Omni-Pure™ Plant DNA Purification System; 100 purifications	40-4060-01	100 Purifications
Omni-Pure™ Genomic Blood DNA Purification System; 100 purifications	40-4010-01	100 Purifications
Omni-Pure™ Genomic Blood DNA Purification System; 500 purifications	40-4010-05	500 Purifications