



Product Specification Summary

Genetic Tools & Reagents

Catalog Number	40-4026-30
Product Name	G3-DNA Wash solution; 2X concentrate; 2X 150 mL
Size	300 mL
Description	G3-DNA Wash solution; 2X concentrate; 2X 150 mL
Component/Note	G3-DNA Wash solution; 2X concentrate; 2X 150 mL

Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17375>



Product Manual http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-4026_G3_V1.1.pdf



Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL