



# Product Specification Summary

## Omni-Pure™ RNA Purification System

|                |  |
|----------------|--|
| Catalog Number | 40-4091-50   |
| Product Name   | Omni-Pure Universal RNA Purification System; 500 purifications |
| Size           | 500 Purifications  |
| Description    | Omni-Pure™ Universal RNA Purification System                   |

Gene Link provides a rapid purification system for extraction of viral DNA or RNA from human bodily fluids including blood. Viral DNA or RNA is captured on a special membrane and then eluted in a low volume for direct use in qualitative and quantitative amplification protocols for detection of a pathogen. The Viral RNA purification system is ideal for small volumes of human bodily fluid samples, i.e., serum, plasma and CSF. Using the easy spin column format, purification of HIV, HCV and other RNA viruses is easily accomplished in less than 30 minutes and ready for RT-PCR amplification.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1084>



Product Manual <http://www.genelink.com/Literature/ps/M40-4081-XX V4.1 RNA Fluid&Spin.pdf>



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_40-Omni-Pure.pdf](http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf)



### Related Products

| Product  | Catalog No | Size              |
|--|------------|-------------------|
| Omni-Pure Viral RNA Spin Column Purification System; 100 purifications | 40-3650-01 | 100 Purifications |
| Omni-Pure Viral RNA Spin Column Purification System 500 purifications  | 40-3650-05 | 500 Purifications |

Omni-Pure Blood & Bodily fluids RNA Purification System; 50 purifications  
190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS40-4091-50.pdf Print Date Version : October 20, 2018 Page 1 of 1

