

Product Specification Summary

Omni-Matrix Sample Collection System

Catalog Number 40-4100-05
Product Name Omni-Mag™ DNA Purification System Kit; 500 purifications.
Size 500 Purifications
Description Omni-Mag DNA Purification System Kit; 500 purifications.

Omni-Mag™ DNA Purification System Kit An optimized kit for DNA extractions from trace samples and specifically touch DNA samples for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17368>



Product Manual http://www.genelink.com/Literature/ps/M40-4100_Omni-Mag_V1.4.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-4100_Omni-Mag.pdf



Related Products

Product	Catalog No	Size
G3 DNA Wash buffer 2X concentrate; 150 mL	40-4026-15	150 mL
Omni-Mag™ Lysis Buffer; 75 mL	40-4119-75	75 mL
Omni-Mag™ DNA Purification System Kit, 100 purifications	40-4100-01	100 Purifications