



Product Specification Summary

Omni-Clean™ DNA Concentration Systems

Catalog Number	40-4112-15
Product Name	Glass Beads; 1.5 mL; Omni-Clean DNA Purification System
Size	1.5 mL
Description	Glass Beads; 1.5 mL; Omni-Clean DNA Purification System
Component/Note	Glass Beads; 1.5 mL; Omni-Clean DNA Purification System

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17247>



Product Manual <http://www.genelink.com/Literature/ps/M40-41XX-XX V3.2.pdf>



Product MSDS http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf



Related Products

Product	Catalog No	Size
Omni-Clean™ DNA Concentration System; 100 purifications	40-4130-10	100 Purifications
Omni-Clean™ DNA Concentration System; 500 purifications	40-4130-50	500 Purifications
Omni-Clean™ DNA Spin Column Concentration System; 100 purification	40-4140-10	100 Purifications