



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number 40-4118-40  
Product Name Omni-Mag™ PMP (Paramagnetic particles); 4 mL  
Size 4 mL  
Description Omni-Mag™ PMP (Paramagnetic particles); 4 mL  
Component/Note Omni-Mag™ PMP (Paramagnetic particles); 4 mL  
Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17373>



Product Manual [http://www.genelink.com/Literature/ps/M40-4100\\_V1.2\\_Omni-Mag.pdf](http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/SDS\\_40-4118\\_PMP\\_V1.2.pdf](http://www.genelink.com/Literature/ps/SDS_40-4118_PMP_V1.2.pdf)



### Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL