



## Product Specification Summary

### Omni-Clean™ Gel DNA Purification Systems

Catalog Number	40-4119-75
Product Name	Omni-Mag™ Lysis Buffer; 75 mL
Size	75 mL
Description	Omni-Mag™ Lysis Buffer; 75 mL
Component/Note	Omni-Mag™ Lysis Buffer; 75 mL

Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17357>



Product Manual [http://www.genelink.com/Literature/ps/M40-4100\\_V1.2\\_Omni-Mag.pdf](http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/SDS\\_40-4119\\_Lysis\\_Buffer\\_V1.2.pdf](http://www.genelink.com/Literature/ps/SDS_40-4119_Lysis_Buffer_V1.2.pdf)



#### Related Products

Product	Catalog No	Size
Omni-Clean™ Gel DNA Beads Purification System; 100 purification	40-4110-10	100 Purifications
Omni-Clean™ Gel DNA Beads Purification System; 500 purifications	40-4110-50	500 Purifications
Omni-Clean™ Gel DNA Spin Column Purification System; 100 purifications	40-4120-10	100 Purifications