



## Product Specification Summary

### Omni-Clean™ DNA Concentration Systems

Catalog Number	40-4130-10
Product Name	Omni-Clean™ DNA Concentration System; 100 purifications
Size	100 Purifications
Description	Omni-Clean™ DNA Concentration System

The Omni-Clean™ DNA Concentration System yields higher quality DNA than with conventional ethanol precipitation. The samples are completely desalted, purified and concentrated. It is ideal for sequencing and genotyping samples with stringent purity requirements. The electropherogram shows samples purified using the Omni-Clean™ system.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=984>



Product Manual <http://www.genelink.com/Literature/ps/M40-41XX-XX V3.2.pdf>



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_40-Omni-Pure.pdf](http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf)



#### Related Products

Product	Catalog No	Size
Omni-Clean™ DNA Concentration System; 100 purifications	40-4130-10	100 Purifications
Omni-Clean™ DNA Concentration System; 500 purifications	40-4130-50	500 Purifications
Omni-Clean™ DNA Spin Column Concentration System; 100 purification	40-4140-10	100 Purifications