



Product Specification Summary

PCR Amplification & Analysis

Catalog Number	40-4140-10
Product Name	Omni-Clean™ DNA Spin Column Concentration System; 100 purification
Size	100 Purifications
Description	Omni-Clean™ DNA Spin Column Concentration System

The Omni-Clean™ DNA Concentration System yields higher quality DNA than with conventional ethanol precipitation. The samples are completely desalted, purified and concentrated. It is ideal for sequencing and genotyping samples with stringent purity requirements. The electropherogram shows samples purified using the Omni-Clean™ system.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=986>



Product Manual <http://www.genelink.com/Literature/ps/M40-41XX-XX V3.2.pdf>



Product MSDS http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf



Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL