



## Product Specification Summary

### Omni-Clean™ DNA Concentration Systems

|                |  |
|----------------|--|
| Catalog Number | 40-4140-50   |
| Product Name   | Omni-Clean™ DNA Spin Column Concentration System; 500 purification |
| Size           | 500 Purifications  |
| Description    | Omni-Clean™ DNA Spin Column Concentration System                   |

The Omni-Clean™ DNA Concentration System yields higher quality DNA than with conventional ethanol precipitation. The samples are completely desalted, purified and concentrated. It is ideal for sequencing and genotyping samples with stringent purity requirements. The electropherogram shows samples purified using the Omni-Clean™ system.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=987>



Product Manual <http://www.genelink.com/Literature/ps/M40-41XX-XX V3.2.pdf>



Product MSDS [http://www.genelink.com/Literature/ps/MSDS\\_40-Omni-Pure.pdf](http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf)



#### Related Products

| Product  | Catalog No | Size              |
|--|------------|-------------------|
| Omni-Clean™ DNA Concentration System; 100 purifications            | 40-4130-10 | 100 Purifications |
| Omni-Clean™ DNA Concentration System; 500 purifications            | 40-4130-50 | 500 Purifications |
| Omni-Clean™ DNA Spin Column Concentration System; 100 purification | 40-4140-10 | 100 Purifications |