

# Product Specification Summary

## Omni-Matrix Sample Collection System

Catalog Number	40-5016-15
Product Name	DNA Elution Buffer (Low TE); 10 X 1.5 mL
Size	15 mL
Description	DNA Elution Buffer (Low TE); 10 X 1.5 mL
Component/Note	DNA Elution Buffer (Low TE); 10 X 1.5 mL

Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17376>



Product Manual [http://www.genelink.com/Literature/ps/M40-4100\\_V1.2\\_Omni-Mag.pdf](http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/SDS\\_40-5016\\_Elution\\_Buffer\\_V1.1.pdf](http://www.genelink.com/Literature/ps/SDS_40-5016_Elution_Buffer_V1.1.pdf)



### Related Products

Product	Catalog No	Size
G3 DNA Wash buffer 2X concentrate; 150 mL	40-4026-15	150 mL
Omni-Mag™ Lysis Buffer; 75 mL	40-4119-75	75 mL
Omni-Mag™ DNA Purification System Kit, 100 purifications	40-4100-01	100 Purifications