



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-5016-15
Product Name	DNA Elution Buffer (Low TE); 10 X 1.5 mL
Size	15 mL
Description	DNA Elution Buffer (Low TE); 10 X 1.5 mL
Component/Note	DNA Elution Buffer (Low TE); 10 X 1.5 mL

Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17376>



Product Manual [http://www.genelink.com/Literature/ps/M40-4100\\_V1.2\\_Omni-Mag.pdf](http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf)



Product MSDS [http://www.genelink.com/Literature/ps/SDS\\_40-5016\\_Eluion\\_Buffer\\_V1.1.pdf](http://www.genelink.com/Literature/ps/SDS_40-5016_Eluion_Buffer_V1.1.pdf)



### Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL