

Product Specification Summary

Omni-Matrix Sample Collection System

Catalog Number	40-5016-30
Product Name	DNA Elution Buffer (Low TE); 30 X 1.5 mL
Size	30 mL
Description	DNA Elution Buffer (Low TE); 30 X 1.5 mL
Component/Note	DNA Elution Buffer (Low TE); 30 X 1.5 mL

Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17377>



Product Manual http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-5016_Elution_Buffer_V1.1.pdf



Related Products

Product	Catalog No	Size
G3 DNA Wash buffer 2X concentrate; 150 mL	40-4026-15	150 mL
Omni-Mag™ Lysis Buffer; 75 mL	40-4119-75	75 mL
Omni-Mag™ DNA Purification System Kit, 100 purifications	40-4100-01	100 Purifications