



Product Specification Summary

Genetic Tools & Reagents

| | |
|----------------|------------------------------------------|
| Catalog Number | 40-5016-30 |
| Product Name | DNA Elution Buffer (Low TE); 30 X 1.5 mL |
| Size | 30 mL |
| Description | DNA Elution Buffer (Low TE); 30 X 1.5 mL |
| Component/Note | DNA Elution Buffer (Low TE); 30 X 1.5 mL |

Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17377>



Product Manual http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-5016_Elution_Buffer_V1.1.pdf



Related Products

| Product | Catalog No | Size |
|----------------------------------------------------|------------|--------|
| Omni-Array™ mRN Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |