



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-5020-01
Product Name	Hybwash A, Hyb Wash Solution Conc. (20X SSC); 1L
Size	1 L
Description	Hybwash A, Hyb Wash Solution Conc. (20X SSC); 1L
Component/Note	Hybwash A, Hyb Wash Solution Conc. (20X SSC); 1L
Storage	Store at 15-25C

Hybwash A, Hyb Wash Solution Conc. (20X SSC); 1L

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17322>



Product Manual [http://www.genelink.com/Literature/ps/Dig\\_Southern\\_Reagents\\_V3.4.pdf](http://www.genelink.com/Literature/ps/Dig_Southern_Reagents_V3.4.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL