



# Product Specification Summary

## Genetic Tools & Reagents

|  |  |
|--|--|
| Catalog Number   | 40-5021-10   |
| Product Name   | Hybwash B, Hyb Wash Solution Concentrate (10% SDS); 100 mL   |
| Size   | 100 mL   |
| Description  | Hybwash B, Hybridization Wash Solution Concentrate (10% SDS) |
| Component/Note   | Hybwash B, Hybridization Wash Solution Concentrate (10% SDS) |
| Component/Note   | Hybwash B, Hybridization Wash Solution Concentrate (10% SDS) |
| Hybwash B, Hybridization Wash Solution Concentrate (10% SDS) |  |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=84>



Product Manual [http://www.genelink.com/Literature/ps/Dig\\_Southern\\_Reagents\\_V3.4.pdf](http://www.genelink.com/Literature/ps/Dig_Southern_Reagents_V3.4.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product   | Catalog No | Size   |
|---|------------|--------|
| Omni-Array™ mRNA Antisense strand amplification kit | 08-0021-02 | 1 Kit  |
| Omni-Array™ Sense Strand mRNA Amplification Kit     | 08-0011-02 | 1 Kit  |
| Omni-Marker™ Universal unlabeled; 100 uL            | 40-3005-01 | 100 uL |