



Product Specification Summary

Genetic Tools & Reagents

| | |
|----------------|--|
| Catalog Number | 40-5023-20 |
| Product Name | Lumisol II, Hybridization Solution; 200 mL |
| Size | 200 mL |
| Component/Note | Ready-to-use hybridization solutions specifically formulated for non-radioactive chemiluminescent hybridization. |

Lumisol™ II hybridization solution can be used for all stringent types of nucleic acid blot hybridization conditions, particularly using non-radioactive digoxigenin labeled probes. Lumisol™ II hybridization solution can be used for prehybridization and hybridization. The hybridization temperature should be calculated for proper results. Duration of hybridization can be reduced to 6 hours and overnight hybridization can be used for high sensitivity requirements and for convenience.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=94>



Product Manual http://www.genelink.com/Literature/ps/PS40-5023_Lumisol_Ver_1.5.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|--------|
| Omni-Array™ mRN Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |