



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-5026-10
Product Name	Blocking solution 10%; 100 mL
Size	100 mL
Description	Blocking solution 10% for hybridization; 100 mL
Component/Note	10% Blocking solution for hybridization
Component/Note	

Blocking solution for hybridization (10%) is a specially formulated reagent containing casein for blocking of Southern blot membranes for chemi-luminescent detection protocols.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1>



Product Manual [http://www.genelink.com/Literature/ps/Dig\\_Southern\\_Reagents\\_V3.4.pdf](http://www.genelink.com/Literature/ps/Dig_Southern_Reagents_V3.4.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	