



Product Specification Summary

Genetic Tools & Reagents

Catalog Number	40-5027-00
Product Name	Seq. Loading buffer; 1 mL
Size	1 mL
Description	2 X Seq Loading Buffer
Component/Note	CautionContains Formamide
Component/Note	DENATURING
Storage	Store at 15-25°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=150>



Product Manual <http://www.genelink.com/Literature/ps/PS40-5027-XXver2.2.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit

Omni-Marker™ Universal unlabeled; 100 uL

40-3005-01

190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |

SS40-5027-00.pdf

Print Date Version : February 23, 2018

Page 1 of 1

