



Product Specification Summary

Genetic Tools & Reagents

| | |
|----------------|--|
| Catalog Number | 40-5027-10 |
| Product Name | Loading Buffer 2X BPB/XC Denaturing for Sequencing; 1 mL |
| Size | 1 mL |
| Description | 2X Concentrate |
| Component/Note | Denaturing for Sequencing |

Gel Loading Buffer 2X BPB/XC denaturing for sequencing is a formulation that contains bromophenol blue (BPB) and xylene cyanol (XC) in formamide. This is specifically used for DNA sequencing sample tracking dyes in denaturing polyacrylamide gels.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=923>



Product Manual <http://www.genelink.com/Literature/ps/PS40-5027-XXver2.2.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|--------|
| Omni-Array™ mRNA Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |