



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-5027-15
Product Name	Loading buffer 2X BPB/XC Sequencing; 15 mL
Size	15 mL
Description	Loading buffer 2X BPB/XC Sequencing
Component/Note	Gel Loading Buffer 2X BPB/XC denaturing for sequencing is a formulation that contains bromophenol blue (BPB) and xylene cyanol (XC) in formamide.
Component/Note	Gel Loading Buffer 2X BPB/XC denaturing for sequencing is a formulation that contains bromophenol blue (BPB) and xylene cyanol (XC) in formamide. This is specifically used for DNA sequencing sample tracking dyes in denaturing polyacrylamide gels.
Component/Note	Gel Loading Buffer 2X BPB/XC denaturing for sequencing is a formulation that contains bromophenol blue (BPB) and xylene cyanol (XC) in formamide. This is specifically used for DNA sequencing sample tracking dyes in denaturing polyacrylamide gels.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=922>



Product Manual <http://www.genelink.com/Literature/ps/PS40-5027-XXver2.2.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit