



Product Specification Summary

Genetic Tools & Reagents

| | |
|---|---|
| Catalog Number | 40-5034-10 |
| Product Name | Depurination Solution (2X) for Southern Blotting; 1 L |
| Size | 1 L |
| Description | Depurination Solution 2X for Southern Blotting;1L |
| Component/Note | Depurination Solution 2X for Southern Blotting;1L |
| Component/Note | Depurination Solution 2X for Southern Blotting;1L |
| Depurination Solution 2X for Southern Blotting;1L | |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17098>



Product Manual http://www.genelink.com/Literature/ps/Dig_Southern_Reagents_V3.4.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|--------|
| Omni-Array™ mRNA Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |