



# Product Specification Summary

## Genetic Tools & Reagents

|   |   |
|---|---|
| Catalog Number                                    | 40-5035-10  |
| Product Name                                      | Denaturation Solution (2X) for Southern Blotting; 1 L |
| Size  | 1 L   |
| Description                                       | Denaturation Solution 2X for Southern Blotting;1L     |
| Component/Note                                    | Denaturation Solution 2X for Southern Blotting;1L     |
| Component/Note                                    | Denaturation Solution 2X for Southern Blotting;1L     |
| Denaturation Solution 2X for Southern Blotting;1L |   |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17099>



Product Manual [http://www.genelink.com/Literature/ps/Dig\\_Southern\\_Reagents\\_V3.4.pdf](http://www.genelink.com/Literature/ps/Dig_Southern_Reagents_V3.4.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product   | Catalog No | Size   |
|---|------------|--------|
| Omni-Array™ mRNA Antisense strand amplification kit | 08-0021-02 | 1 Kit  |
| Omni-Array™ Sense Strand mRNA Amplification Kit     | 08-0011-02 | 1 Kit  |
| Omni-Marker™ Universal unlabeled; 100 uL            | 40-3005-01 | 100 uL |