

Product Specification Summary

DNA & RNA Reconstitution Solutions

Catalog Number	40-5106-15
Product Name	DNA Elution Buffer (Low TE); 10 X 1.5 mL
Size	15 mL
Description	DNA Elution Buffer (Low TE)
Component/Note	DNA Elution Buffer (Low TE) solution consists of 10 mM Tris pH 8.0 and 0.1 mM EDTA.

DNA Elution Buffer (Low TE) solution consists of 10 mM Tris pH 8.0 and 0.1 mM EDTA. This is sterile and filtered through a 0.2 µm filter. All Gene Link TE Buffers are provided nuclease-free, ready for use, and require no extra preparation. Alkalinity and free cations have been implicated in catalyzing the hydrolysis of RNA. TE Buffer pH 7.0 reduces base hydrolysis through chelation of divalent cations by EDTA and through resistance to pH changes by the Tris buffer.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17369>



Product Manual http://www.genelink.com/Literature/ps/PS40-30-DNA&RNA_Soln_Ver2.1.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDS_40-DNA&RNA_solutions.pdf



Related Products

Product	Catalog No	Size
Nuclease Free Water (DEPC Free) 10 X 1.6 mL; 16 mL	40-3001-16	16 mL
DEPC Treated Water; 10 X 1.6 mL	40-3000-16	1.6 mL