



# Product Specification Summary

## PCR Amplification & Analysis

|                |   |
|----------------|---|
| Catalog Number | 40-5200-40  |
| Product Name   | Taq DNA Polymerase 400 units 5 u/ul; 80 uL; 400 units |
| Size           | 400 each  |
| Description    | Taq DNA Polymerase 400 units 5 u/ul; 80 uL; 400 units |
| Component/Note | 5 u/uL  |
| Storage        | Store at -20°C  |

Taq DNA Polymerase 400 units 5 u/ul; 80 uL; 400 units

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1053>



Product Manual [http://www.genelink.com/Literature/ps/M40-5200-XX\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/M40-5200-XX_Ver3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product                                  | Catalog No | Size   |
|--|------------|--------|
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |
| Omni-Marker™ Universal unlabeled; 1 mL   | 40-3005-10 | 1 mL   |
| Omni-Marker™ Low unlabeled; 500 uL       | 40-3006-05 | 500 uL |