



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-5200-40
Product Name	Taq DNA Polymerase 400 units 5 u/ul; 80 uL; 400 units
Size	400 each
Description	Taq DNA Polymerase 400 units 5 u/ul; 80 uL; 400 units
Component/Note	5 u/uL
Storage	Store at -20°C

Taq DNA Polymerase 400 units 5 u/ul; 80 uL; 400 units

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1053>



Product Manual [http://www.genelink.com/Literature/ps/M40-5200-XX\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/M40-5200-XX_Ver3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL