



Product Specification Summary


Omni-Pure™ Genomic DNA Purification Systems

| | |
|----------------|---|
| Catalog Number | 40-5205-01 |
| Product Name | Lytic Enzyme 4000 units/mL; 1 mL |
| Size | 1 mL |
| Description | Lytic Enzyme 4000 units/mL; 1 mL |
| Component/Note | Component of OmniPure DNA Purification System Kit |
| Component/Note | From <i>Arthrobacter luteus</i> |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17496>

Product Manual  http://www.genelink.com/Literature/ps/GDComponents_40-4010-XX.pdf

Product MSDS  http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf

Related Products

| Product | Catalog No | Size |
|---|------------|-------------------|
| Omni-Pure™ Plant DNA Purification System; 100 purifications | 40-4060-01 | 100 Purifications |
| Omni-Pure™ Genomic Blood DNA Purification System; 100 purifications | 40-4010-01 | 100 Purifications |
| Omni-Pure™ Genomic Blood DNA Purification System; 500 purifications | 40-4010-05 | 500 Purifications |