



Product Specification Summary

PCR Amplification & Analysis

Catalog Number 40-5211-01
Product Name Taq PCR Kit; 200 rxns
Size 200 rxns
Description PCR Products

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1054>



Product Manual http://www.genelink.com/Literature/ps/M40-5200-XX_Ver3.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|--------|
| Omni-Marker™ Universal unlabeled; 100 uL | 40-3005-01 | 100 uL |
| Omni-Marker™ Universal unlabeled; 1 mL | 40-3005-10 | 1 mL |
| Omni-Marker™ Low unlabeled; 500 uL | 40-3006-05 | 500 uL |